

Contents:

- 1 aluminium container with 100 test sticks
- 2 bottles Molybdenum-1
- 1 measuring tube with 5 ml mark
- 1 measuring spoon

Hazard warning:

This test does not contain any amounts of harmful substances which must be specially labelled as hazardous.

Method of application:

1. Rinse measuring tube with the sample and fill to the 5 ml mark.
2. Add **2 level measuring spoons** of Molybdenum-1 and dissolve by swirling carefully.
3. Remove only as many test sticks as are required, and reseal container immediately after use. Do not touch the test paper zone.
4. Place test stick for 1 min. into prepared sample solution.
5. Remove test stick, shake off excess liquid and compare to colour scale.

Remark:

The pH value of the test solution should be below pH 1 after adding Molybdenum-1. If the pH value is not reached with the quantity of the reagent Molybdenum-1 adjust the pH value to 4-6 by adding further reagent Molybdenum-1, or, prior to the determination with diluted sulfuric acid.

Storage:

Avoid exposing the sticks to sunlight and moisture. Store the kit below +30 °C in a dry place.