visocolor<sup>®</sup> alpha

Ammonio  
Amonio

Ammonium

0.2 - 3 mg/l NH<sub>4</sub><sup>+</sup>

Art.-Nr. 935 012

50 Tests

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on the chlorine may reduce the measurement reading or suppress the reaction totally. ions and produce higher results. Depending on their concentration, substances which draw test using sodium hydroxide solution. Primary amines react in the same way as ammonium

**Interférences:** Strongly acidic or buffered solutions have to be adjusted to pH 9-10 for the test using sodium hydroxide solution. Primary amines react in the same way as ammonium ions and produce higher results. Depending on their concentration, substances which draw test using sodium hydroxide solution. Primary amines react in the same way as ammonium

**Dispensing of the sample:** The used analysis specimens can be flushed down the drain with tap water and channelled off to the local sewage treatment works. This method can be used for analysing sea water after dilution (1:10). After 5 min place the measuring vessel on the colour chart and assign the value by comparison of the colour. Mid-values can be estimated.

**Instructions for use:** also refer to the pictogram on the inside. This test does not contain any amounts of harmful substances which must be specially labelled as hazardous.

**Test kit for performing colorimetric tests on ammonium**

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